

Product Information Specifications & Tests

DL-alpha-Tocopherol

Product code: 0410276

1. **Appearance:** **clear, viscous oil**
Visual evaluation

2. **Colour** **colourless to yellowish-brown**
Visual evaluation

3. **Identity** **corresponds**
 - Identification Test A. (Optical rotation) as described in the Ph. Eur. monograph 0692 of all-*rac*- α -Tocopherol.
 - Gas chromatography method as described under **Assay (Ph. Eur.)**: Identity corresponds when the retention time of the main peak from the sample chromatogram matches the retention time of main peak from the reference solution.

4. **Optical Rotation** **-0.01° to + 0.01°**
Proceed according to Ph. Eur. method Optical Rotation (2.2.7) and Identification test A. as described in the Ph. Eur. Monograph 0692 of all-*rac*- α -Tocopherol.

5. **Refractive Index 589 nm, 20°C** **1.503 - 1.507**
Proceed according to Ph. Eur. method Refractive Index (2.2.6).

6. **Absorbance in Ethanol at about 292 nm (Max.)** **71.0 - 76.0**
Proceed according to Ph. Eur. method Absorption Spectrophotometry, Ultraviolet and Visible (2.2.25).

Stauber Item # 20790

Product Information

Specifications & Tests

DL-alpha-Tocopherol

Product code: 0410276

7. Acidity

max. 1.0 mL 0.10 N NaOH

Proceed according to the specific test for Acidity as described in the USP monograph of Vitamin E or with the in-house as described, below:

Reagents, solvents, standard materials

Ethanol, abs.	96 % Ph. Eur. <i>R</i>
Diethyl ether	≥99.7 % Ph. Eur. <i>R</i>
1 N NaOH	Ph. Eur. <i>R</i>

The described system or a system with equivalent performance may be used:

Titration System Titration unit for potentiometric titration, e.g., 904 Titrand (Metrohm) with 5 mL dosing unit and stirrer, or equiv.

Electrode: pH electrode for non-aqueous acid/base titrations, e.g., Solvotrode (Metrohm), or equiv.

Analysis of Sample Dissolve 2-3 g of the sample in 50 mL of a mixture of equal volumes of ethanol (96 per cent) *R* and diethyl ether *R*. Titrate while stirring with 0.02 mol/L NaOH, using the Solvotrode. If necessary, heat to about 40 °C to dissolve the sample. When heating has been applied to aid dissolution, maintain the temperature at about 40 °C during the titration.

Analysis of Blank Sample Proceed as described above with a blank sample containing 50 mL of a mixture of equal volumes of ethanol (96 per cent) *R* and diethyl ether *R*.

Calculation

$$\text{Acidity} = \frac{(M - B) \times 0.2}{E}$$

M:	Consumption of 0.02 N NaOH for the sample solution [mL]
B:	Consumption of 0.02 N NaOH for the blank solution [mL]
0.2:	Conversion from 0.02 N NaOH to 0.1 N NaOH
E:	Sample Weight [g]

Product Information

Specifications & Tests

DL-alpha-Tocopherol

Product code: 0410276

8. Sulfated Ash (Residue on Ignition) max. 0.1 %

Proceed according to Ph. Eur. method Sulfated Ash (2.4.14).

9. Residual Solvents

Toluene max. 80 ppm

Determine according to Ph. Eur. method Identification and control of residual solvents (2.4.24) or with in-house gas chromatographic method as described, below:

Reagents, solvents, standard materials

DMAC N,N-Dimethylacetamide ≥ 99.9 %

Toluene ≥ 99.7 %

Preparation of test solution, blank solution, and calibration solution

Stock solution In a 50 mL volumetric flask, add approx. 20 mL DMAC, then add 890 mg toluene. Fill up to volume with dimethylacetamide.

Blank solution DMAC

SST solution In a 50 mL volumetric flask, add approx. 20 mL DMAC, then add 250 μ L stock solution. Fill up to volume with dimethylacetamide and homogenize. Transfer 2 mL into a 20 mL headspace vial and seal hermetically.

Test solution Transfer approx. 200 mg of the test substance, accurately weighted, into a 20 mL Headspace-Vial, add 2 mL of DMAC and seal hermetically.

SST Requirements

Blank No interfering peak at toluene retention time.

SST The user laboratory should have procedures in place for adequate system suitability testing.

Product Information Specifications & Tests

DL-alpha-Tocopherol

Product code: 0410276

Equipment

The described system or a system with equivalent performance may be used

GC system	gas chromatography with static head-space injection
Stationary phase	Agilent DB-1 GC (dimethylpolysiloxane): length 60 m, ID 0.32 mm, Film thickness 3.0 µm or equivalent

GC & HS parameters

Static head space (HS) sample preparation	Vial equilibration temperature	90 °C (HS oven temperature)	
	Vial equilibration time	15 min	
	Vial Pressurization equilibration time	0.1 min	
	Vial shaking	Low	
	Vial pressure	1.5 bar	
	Loop temperature	110 °C	
	Transfer line temperature	130 °C	
Other HS parameters (for information)	Vial pressurization time	0.02 min	
	Loop fill time	0.05 min	
	Loop fill equilibration time	0.05 min	
	Temperature equilibration time	0.2 min	
	Sampler inject time	2 min	
	Sequence valve purge time	2 min	
Oven temperature	Rate (°C/min)	Temperature (°C)	Isotime (min)
	---	60	5
	30	120	6
	22.5	210	2

Product Information Specifications & Tests

DL-alpha-Tocopherol

Product code: 0410276

GC parameters	Carrier gas	Helium for chromatography
	Carrier flow	Constant flow, 1.5 ml/min
	Split	1:10
	Injection volume	1 ml
	Injector temperature	150 °C
	Detector temperature	250 °C
	Chromatogram run time	19 min

Calculation

According to the determination by area using an external standard

$$\text{Content [ppm]} = \frac{A_S \times W_{\text{Std}} \times CF}{A_{\text{Std}} \times W_S} \times 200$$

A_S = area of the corresponding residual solvent in the sample solution [pA*min]

W_{Std} = amount of the corresponding residual solvent in the reference solution [mg]

A_{Std} = area of the corresponding residual solvent in the reference solution [pA*min]

W_S = amount of the test substance [mg]

CF = correction factor

200 = calculation factor in ppm (10^6) with dilution factor (2×10^{-4})

Retention time (for information)

Toluene: 15.3 min

Limit of quantification (LOQ): 20 ppm

Product Information Specifications & Tests

DL-alpha-Tocopherol

Product code: 0410276

10. Related Substances (Ph. Eur.)

• Impurity A	max. 0.5 %
• Impurity B	max. 1.5 %
• Impurities C and D	max. 1.0 %
• Any other impurity	each max. 0.25 %
• Total	max. 2.5 %

Proceed according to test “Related substances” as described in the Ph. Eur. monograph 0692 of all-*rac*- α -Tocopherol or to the in-house gas chromatographic method as described under Assay (Ph. Eur.).

11. Assay (Ph. Eur.)

97.0 - 102.0 %

Proceed according to test “Assay” as described in the Ph. Eur. monograph 0692 of all-*rac*- α -Tocopherol or the the gas chromatographic in-house method as described, below.

The described system or a system with equivalent performance may be used:

Gas chromatograph	Agilent, Model 6850 with FID detector
Column	Fused silica (crosslinked 5 % Phenylsilanol, 95 % Methylsilanol)
Dimensions:	25 m x 0.32 mm
Film thickness:	0.52 μ m
Carrier gas:	H ₂
Total flow rate:	2.5 mL/min.
Split ratio:	1:5
Column temperature:	Start 180 °C (2 min.), increase 20.0 °C/min. to 280 °C (12 min.)
Injector temperature:	280 °C
Detector temperature:	300 °C
Injection volume:	1.0 μ L
Duration:	19 min.

Product Information

Specifications & Tests

DL-alpha-Tocopherol

Product code: 0410276

Solutions

Internal standard solution

Weigh in 33.9 to 34.1 g C₂₈ alkane (= internal standard) and dissolve with 10 L heptane.

Sample solution

Weigh in 170.0 to 210.0 mg of sample (accurately weighed) into a 50 mL volumetric flask and make up to volume with the internal standard solution. Mix for 4 min. (injection solution).

Reference solution

Weigh in 170.0 to 210.0 mg (accurately weighed) of the current all-*rac*-α-tocopherol reference standard into a 50 mL volumetric flask, make up to volume with internal standard solution. Mix for 4 min. (injection solution).

Identity

Identity corresponds when the retention time of the main peak from the sample chromatogram matches the retention time of main peak from the reference solution.

Retention time [min.]	Compound	Relative retention time
approx. 9.0	Internal standard	0.71
approx. 10.5	Impurity A	0.83
approx. 11.4	Impurity B	0.90
approx. 12.7	all- <i>rac</i> -α-tocopherol	1.00
approx. 13.0	Impurity C + D	1.02

Calculation

For calculating the amount of impurity A, B, C, and D, the response factor of all-*rac*-α-tocopherol is used.

Product Information Specifications & Tests

DL-alpha-Tocopherol

Product code: 0410276

Calculation of the assay:

$$\text{Assay (\%, w/w)} = \frac{A_S \times RF \times W_{iStd(S)} \times 100}{A_{iStd(S)} \times W_S}$$

With RF:

$$RF = \frac{W_{Ref} \times A_{Ph.Eur.} \times A_{iStd(St)}}{W_{iStd(St)} \times A_{St} \times 100\%}$$

A_S	peak area of all- <i>rac</i> - α -tocopherol in sample solution chromatogram
A_{St}	peak area of all- <i>rac</i> - α -tocopherol in standard solution chromatogram
$A_{iStd(S)}$	peak area of internal standard (Octacosan) in sample chromatogram
$A_{iStd(St)}$	peak area of internal standard (Octacosan) in standard solution chromatogram
$A_{Ph.Eur.}$	assay of Ph. Eur. reference standard [%]
$W_{iStd(S)}$	amount of internal standard in sample solution [mg]
$W_{iStd(St)}$	amount of internal standard in standard solution [mg]
W_S	sample weight [mg]
W_{Ref}	weight of all- <i>rac</i> - α -tocopherol reference standard [mg]
RF	response factor all- <i>rac</i> - α -tocopherol
100%	calculation factor to percent

This document is valid without signature.

Disclaimer: The information contained herein is to the best of our knowledge accurate. However, we refuse any liability for the application and use of further processed material containing our product. Solely the producer of the final product has to assume full responsibility according to the relevant local regulations. The content of this document is subject to change without further notice unless otherwise agreed in writing. Please consult your local representative or one of DSM's websites at <https://webshop.dsm.com> for the latest version of this document or the latest information.